

NCRI Lymphoma & Radiotherapy Groups QA Programme

TLDs FOR ELECTRONS

Instructions:

Positioning of TLD holders for patient measurements.

TLDs should be placed at the centre of the treatment field. If using bolus the TLDs should be placed under the bolus. The TLDs should be secured in position prior to treatment and not removed until after treatment for that day.

Complete the following information.

Treatment area	
Dose per fraction	
Number of fractions	
Electron energy	
Field size	
Expected skin surface dose	
Anticipated 90% depth	
Thickness of bolus	

Please attach a hard copy of the patients treatment plan if available (central axis only).

Centre	
Trial No	
Patient Initials	
Date of Measurement	
Machine Output for that Day	

Please return the TLDs even if they are not used or if they are used incorrectly.